

CNN2 Antibody, HRP conjugated



PACO48727

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thin filament-associated protein that is implicated in the regulation and modulation of smooth muscle contraction. It is capable of binding to actin, calmodulin, troponin C and tropomyosin. The interaction of calponin with actin inhibits the actomyosin Mg-ATPase activity.

Gene ID:

CNN2

Uniprot

Q99439

Synonyms:

Calponin-2 (Calponin H2, smooth muscle) (Neutral calponin), CNN2

Immunogen:

Recombinant Human Calponin-2 protein (2-309AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A