

USP25 Antibody, HRP conjugated



PACO48763

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Deubiquitinating enzyme that hydrolyzes ubiquitin moieties conjugated to substrates and thus, functions to process newly synthesized Ubiquitin, to recycle ubiquitin molecules or to edit polyubiquitin chains and prevents proteasomal degradation of substrates. Hydrolyzes both 'Lys-48'- and 'Lys-63'-linked tetraubiquitin chains.

Gene ID:

USP25

Uniprot

Q9UHP3

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 25 (EC 3.4.19.12) (Deubiquitinating enzyme 25) (USP on chromosome 21) (Ubiquitin thioesterase 25) (Ubiquitin-specific-processing protease 25), USP25, USP21

Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase 25 protein (348-606AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A