

# SLCO3A1 Antibody, HRP conjugated



PACO48771

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Mediates the Na(+)-independent transport of organic anions such as estrone-3-sulfate. Mediates transport of prostaglandins (PG) E1 and E2, thyroxine (T4), deltorphin II, BQ-123 and vasopressin, but not DPDPE (a derivative of enkephalin lacking an N-terminal tyrosine residue), estrone-3-sulfate, taurocholate, digoxin nor DHEAS.

**Gene ID:**

SLCO3A1

**Uniprot**

Q9UIG8

**Synonyms:**

Solute carrier organic anion transporter family member 3A1 (OATP3A1) (Organic anion transporter polypeptide-related protein 3) (OATP-RP3) (OATPRP3) (Organic anion-transporting polypeptide D) (OATP-D) (PGE1 transporter) (Sodium-independent organic anion transporter D) (Solute carrier family 21 member 11), SLCO3A1, OATP3A1 OATPD SLC21A11

**Immunogen:**

Recombinant Human Solute carrier organic anion transporter family member 3A1 protein (439-539AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A