

AJAP1 Antibody, HRP conjugated



PACO48799

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in cell adhesion and cell migration.

Gene ID:

AJAP1

Uniprot

Q9UKB5

Synonyms:

Adherens junction-associated protein 1 (Membrane protein shrew-1), AJAP1, MOT8 SHREW1

Immunogen:

Recombinant Human Adherens junction-associated protein 1 protein (308-394AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A