

# AJAP1 Antibody, Biotin conjugated



PACO48801

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a role in cell adhesion and cell migration.

**Gene ID:**

AJAP1

**Uniprot**

Q9UKB5

**Synonyms:**

Adherens junction-associated protein 1 (Membrane protein shrew-1), AJAP1, MOT8 SHREW1

**Immunogen:**

Recombinant Human Adherens junction-associated protein 1 protein (308-394AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A