

BTRC Antibody, Biotin conjugated



PACO48817

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

BTRC

Uniprot

Q9Y297

Synonyms:

F-box/WD repeat-containing protein 1A (E3RSIkappaB) (Epididymis tissue protein Li 2a) (F-box and WD repeats protein beta-TrCP) (plkappaBalphalpha-E3 receptor subunit), BTRC, BTRCP FBW1A FBXW1A

Immunogen:

Recombinant Human F-box/WD repeat-containing protein 1A protein (397-448AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A