

ACYP2 Antibody, HRP conjugated



PACO48819

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Its physiological role is not yet clear.

Gene ID:

ACYP2

Uniprot

P14621

Synonyms:

Acylphosphatase-2 (EC 3.6.1.7) (Acylphosphatase, muscle type isozyme) (Acylphosphate phosphohydrolase 2), ACYP2, ACYP

Immunogen:

Recombinant Human Acylphosphatase-2 protein (28-87AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A