

# ALDH6A1 Antibody, Biotin conjugated



PACO48829

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a role in valine and pyrimidine metabolism. Binds fatty acyl-CoA.

**Gene ID:**

ALDH6A1

**Uniprot**

Q02252

**Synonyms:**

Methylmalonate-semialdehyde dehydrogenase [acylating], mitochondrial (MMSDH) (Malonate-semialdehyde dehydrogenase [acylating]) (EC 1.2.1.18) (EC 1.2.1.27) (Aldehyde dehydrogenase family 6 member A1), ALDH6A1, MMSDH

**Immunogen:**

Recombinant Human Methylmalonate-semialdehyde dehydrogenase [acylating], mitochondrial protein (274-399AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A