

# GFRA1 Antibody, HRP conjugated



PACO48923

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for GDNF. Mediates the GDNF-induced autophosphorylation and activation of the RET receptor.

**Gene ID:**

GFRA1

**Uniprot**

P56159

**Synonyms:**

GDNF family receptor alpha-1 (GDNF receptor alpha-1) (GDNFR-alpha-1) (GFR-alpha-1) (RET ligand 1) (TGF-beta-related neurotrophic factor receptor 1), GFRA1, GDNFRA, RETL1, TRNR1

**Immunogen:**

Recombinant Human GDNF family receptor alpha-1 protein (230-385AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A