

Glp1r Antibody, Biotin conjugated



PACO48937

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

Glp1r

Uniprot

O35659

Synonyms:

Glucagon-like peptide 1 receptor (GLP-1 receptor) (GLP-1-R) (GLP-1R), Glp1r

Immunogen:

Recombinant Mouse Glucagon-like peptide 1 receptor protein (22-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A