

PACO48942

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May constitute a glutathione peroxidase-like protective system against peroxide damage in sperm membrane lipids. Since the purified porcine enzyme has very little activity towards hydrogen peroxide or organic hydroperoxides the protective effect is not likely to be exerted by its enzymatic activity. Instead, may protect sperm from premature acrosome reaction in the epididymis by binding to lipid peroxides, which might otherwise interact with phospholipase A2 and induce the acrosome reaction.

Gene ID:

GPX5

Uniprot

O18994

Synonyms:

Epididymal secretory glutathione peroxidase (EC 1.11.1.9) (Epididymis-specific glutathione peroxidase-like protein) (EGLP) (Glutathione peroxidase 5) (GPx-5) (GSHPx-5), GPX5

Immunogen:

Recombinant Pig Epididymal secretory glutathione peroxidase protein (22-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A