

KPNA4 Antibody, FITC conjugated



PACO48996

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KPNA4

Uniprot

O00629

Synonyms:

Importin subunit alpha-3 (Importin alpha Q1) (Qip1) (Karyopherin subunit alpha-4), KPNA4, QIP1

Immunogen:

Recombinant Human Importin subunit alpha-3 protein (51-239AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A