

# Lum Antibody, HRP conjugated



PACO48999

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular space, proteinaceous extracellular matrix, collagen binding, cartilage development, response to growth factor, response to organic cyclic compound.

**Gene ID:**

Lum

**Uniprot**

P51886

**Synonyms:**

Lumican (Keratan sulfate proteoglycan lumican) (KSPG lumican), Lum, Lcn Ldc

**Immunogen:**

Recombinant Rat Lumican protein (19-338AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A