

MAP4 Antibody, FITC conjugated



PACO49004

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Non-neuronal microtubule-associated protein. Promotes microtubule assembly.

Gene ID:

MAP4

Uniprot

P27816

Synonyms:

Microtubule-associated protein 4 (MAP-4), MAP4

Immunogen:

Recombinant Human Microtubule-associated protein 4 protein (575-818AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A