

MAP4 Antibody, Biotin conjugated



PACO49005

Product Information

Size:	Protein Background:
50ug	Non-neuronal microtubule-associated protein. Promotes microtubule assembly.
Reactivity:	Gene ID:
Human	MAP4
Source:	Uniprot
Rabbit	P27816
Isotype:	Synonyms:
IgG	Microtubule-associated protein 4 (MAP-4), MAP4
Applications:	Immunogen:
ELISA	Recombinant Human Microtubule-associated protein 4 protein (575-818AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A