

PACO49014

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the degradation of several amino acid, , odd-chain fatty acid, and cholesterol via propionyl-CoA to the tricarboxylic acid, cycle.

Gene ID:

Mut

Uniprot

P16332

Synonyms:

Methylmalonyl-CoA mutase, mitochondrial (MCM) (EC 5.4.99.2) (Methylmalonyl-CoA isomerase), Mut

Immunogen:

Recombinant Mouse Methylmalonyl-CoA mutase, mitochondrial protein (31-748AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A