

TECPR2 Antibody



PACO49058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Probably plays a role as positive regulator of autophagy.

Gene ID:

TECPR2

Uniprot

O15040

Synonyms:

Tectonin beta-propeller repeat-containing protein 2 (WD repeat-containing protein KIAA0329/KIAA0297), TECPR2, KIAA0297 KIAA0329


Immunogen:

Recombinant Human Tectonin beta-propeller repeat-containing protein 2 protein (516-783AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO49058 at dilution of 1:100.



Immunofluorescent analysis of PC-3 cells using PACO49058 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).