

VGF Antibody, FITC conjugated



PACO49080

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the regulation of cell-cell interactions or in synaptogenesis during the maturation of the nervous system.

Gene ID:

VGF

Uniprot

O15240

Synonyms:

Neurosecretory protein VGF [Cleaved into: Neuroendocrine regulatory peptide-1 (NERP-1); Neuroendocrine regulatory peptide-2 (NERP-2); Antimicrobial peptide VGF[554-577]], VGF

Immunogen:

Recombinant Human Neurosecretory protein VGF protein (23-206AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A