

# FEZ1 Antibody, Biotin conjugated



PACO49089

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in axonal outgrowth as component of the network of molecules that regulate cellular morphology and axon guidance machinery. Able to restore partial locomotion and axonal fasciculation to *C. elegans* unc-76 mutants in germline transformation experiments. May participate in the transport of mitochondria and other cargos along microtubules.

**Gene ID:**

FEZ1

**Uniprot**

Q99689

**Synonyms:**

Fasciculation and elongation protein zeta-1 (Zygin I) (Zygin-1), FEZ1

**Immunogen:**

Recombinant Human Fasciculation and elongation protein zeta-1 protein (12-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A