

APOE Antibody, Biotin conjugated



PACO49109

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Mediates the binding, internalization, and catabolism of lipoprotein particles. It can serve as a ligand for the LDL (apo B/E) receptor and for the specific apo-E receptor (chylomicron remnant) of hepatic tissues.

Gene ID:

APOE

Uniprot

P02649

Synonyms:

Apolipoprotein E (Apo-E), APOE

Immunogen:

Recombinant Human Apolipoprotein E protein (21-110AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A