

IGFBP7 Antibody, FITC conjugated



PACO49140

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds IGF-I and IGF-II with a relatively low affinity. Stimulates prostacyclin (PGI₂) production. Stimulates cell adhesion.

Gene ID:

IGFBP7

Uniprot

Q16270

Synonyms:

Insulin-like growth factor-binding protein 7 (IBP-7) (IGF-binding protein 7) (IGFBP-7) (IGFBP-rP1) (MAC25 protein) (PGI₂-stimulating factor) (Prostacyclin-stimulating factor) (Tumor-derived adhesion factor) (TAF), IGFBP7, MAC25 PSF

Immunogen:

Recombinant Human Insulin-like growth factor-binding protein 7 protein (64-171AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A