

MBP Antibody



PACO49146

Product Information

Size:

50ug

Reactivity:

Macaca fascicularis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MBP

Uniprot

G7PWA0

Synonyms:

Myelin basic protein, MBP, EGM_08965

Immunogen:

Recombinant Macaca fascicularis Myelin basic protein (1-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A