

# LRP8 Antibody, FITC conjugated



PACO49184

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LRP8

**Uniprot**

Q14114

**Synonyms:**

Low-density lipoprotein receptor-related protein 8 (LRP-8) (Apolipoprotein E receptor 2), LRP8, APOER2

**Immunogen:**

Recombinant Human Low-density lipoprotein receptor-related protein 8 protein (131-332AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A