

LRP8 Antibody, Biotin conjugated



PACO49185

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

LRP8

Uniprot

Q14114

Synonyms:

Low-density lipoprotein receptor-related protein 8 (LRP-8) (Apolipoprotein E receptor 2), LRP8, APOER2

Immunogen:

Recombinant Human Low-density lipoprotein receptor-related protein 8 protein (131-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A