

# ARHGAP11B Antibody, FITC conjugated



PACO49192

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hominin-specific protein that promotes development and evolutionary expansion of the brain neocortex. Able to promote amplification of basal progenitors in the subventricular zone, producing more neurons during fetal corticogenesis. Does not possess GTPase activator activity.

**Gene ID:**

ARHGAP11B

**Uniprot**

Q3KRB8

**Synonyms:**

Rho GTPase-activating protein 11B (Protein FAM7B1) (Rho-type GTPase-activating protein 11B), ARHGAP11B, FAM7B1

**Immunogen:**

Recombinant Human Rho GTPase-activating protein 11B protein (99-231AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A