

ARHGAP11B Antibody, Biotin conjugated



PACO49193

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Hominin-specific protein that promotes development and evolutionary expansion of the brain neocortex. Able to promote amplification of basal progenitors in the subventricular zone, producing more neurons during fetal corticogenesis. Does not possess GTPase activator activity.

Gene ID:

ARHGAP11B

Uniprot

Q3KRB8

Synonyms:

Rho GTPase-activating protein 11B (Protein FAM7B1) (Rho-type GTPase-activating protein 11B), ARHGAP11B, FAM7B1

Immunogen:

Recombinant Human Rho GTPase-activating protein 11B protein (99-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A