

PACO49198

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probably mediates the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events.

Gene ID:

HDT2

Uniprot

Q56WH4

Synonyms:

Histone deacetylase HDT2 (HD-tuins protein 2) (Histone deacetylase 2b), HDT2, HD2 HD2B HDA4

Immunogen:

Recombinant Arabidopsis thaliana Histone deacetylase HDT2 protein (1-306AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A