

SLC16A9 Antibody, Biotin conjugated



PACO49217

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Proton-linked monocarboxylate transporter. May catalyze the transport of monocarboxylates across the plasma membrane.

Gene ID:

SLC16A9

Uniprot

Q7RTY1

Synonyms:

Monocarboxylate transporter 9 (MCT 9) (Solute carrier family 16 member 9), SLC16A9, C10orf36 MCT9

Immunogen:

Recombinant Human Monocarboxylate transporter 9 protein (185-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A