

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Chromatin regulator which modulates histone acetylation and gene expression in neural precursor cells. May recruit histone deacetylases (HDACs) to the p160 coactivators/nuclear receptor complex to inhibit ligand-dependent transactivation. Has a role in proliferation and development of cortical neural precursors. May also regulate bone homeostasis.

Gene ID:

ANKRD11

Uniprot

Q6UB99

Synonyms:

Ankyrin repeat domain-containing protein 11 (Ankyrin repeat-containing cofactor 1), ANKRD11, ANCO1

Immunogen:

Recombinant Human Ankyrin repeat domain-containing protein 11 protein (1002-1276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A