

GPR180 Antibody



PACO49222

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

extracellular exosome.

Gene ID:

GPR180

Uniprot

Q86V85

Synonyms:

Integral membrane protein GPR180 (Intimal thickness-related receptor), GPR180, ITR


Immunogen:

Recombinant Human Integral membrane protein GPR180 protein (114-168AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO49222 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).