

# WDFY1 Antibody, HRP conjugated



PACO49235

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Positively regulates TLR3- and TLR4-mediated signaling pathways by bridging the interaction between TLR3 or TLR4 and TICAM1. Promotes TLR3/4 ligand-induced activation of transcription factors IRF3 and NF-kappa-B, as well as the production of IFN-beta and inflammatory cytokines.

**Gene ID:**

WDFY1

**Uniprot**

Q8IWB7

**Synonyms:**

WD repeat and FYVE domain-containing protein 1 (FENS-1) (Phosphoinositide-binding protein 1) (WD40- and FYVE domain-containing protein 1) (Zinc finger FYVE domain-containing protein 17), WDFY1, KIAA1435 WDF1 ZFYVE17

**Immunogen:**

Recombinant Human WD repeat and FYVE domain-containing protein 1 protein (261-357AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A