

# PNKD Antibody, HRP conjugated



PACO49251

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Probable hydrolase that plays an aggravative role in the development of cardiac hypertrophy via activation of the NF-kappa-B signaling pathway.

**Gene ID:**

PNKD

**Uniprot**

Q8N490

**Synonyms:**

Probable hydrolase PNKD (EC 3) (Myofibrillogenesis regulator 1) (MR-1) (Paroxysmal nonkinesinogenic dyskinesia protein) (Trans-activated by hepatitis C virus core protein 2), PNKD, KIAA1184 MR1 TAHCCP2

**Immunogen:**

Recombinant Human Probable hydrolase PNKD protein (78-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A