

PNKD Antibody, Biotin conjugated



PACO49253

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable hydrolase that plays an aggravative role in the development of cardiac hypertrophy via activation of the NF-kappa-B signaling pathway.

Gene ID:

PNKD

Uniprot

Q8N490

Synonyms:

Probable hydrolase PNKD (EC 3) (Myofibrillogenesis regulator 1) (MR-1) (Paroxysmal nonkinesigenic dyskinesia protein) (Trans-activated by hepatitis C virus core protein 2), PNKD, KIAA1184 MR1 TAHCCP2

Immunogen:

Recombinant Human Probable hydrolase PNKD protein (78-243AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A