

GRHL3 Antibody



PACO49258

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

GRHL3

Uniprot

Q8TE85

Synonyms:

Grainyhead-like protein 3 homolog (Sister of mammalian grainyhead) (Transcription factor CP2-like 4), GRHL3, SOM TFCP2L4


Immunogen:

Recombinant Human Grainyhead-like protein 3 homolog protein (31-182AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of PC-3 cells using PACO49258 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).