

ALG13 Antibody, FITC conjugated



PACO49324

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Isoform 1: Possible multifunctional enzyme with both glycosyltransferase and deubiquitinase activities.

Gene ID:

ALG13

Uniprot

Q9NP73

Synonyms:

Putative bifunctional UDP-N-acetylglucosamine transferase and deubiquitinase ALG13 (EC 2.4.1.141) (EC 3.4.19.12) (Asparagine-linked glycosylation 13 homolog) (Glycosyltransferase 28 domain-containing protein 1) (UDP-N-acetylglucosamine transferase subunit ALG13 homolog), ALG13, CXorf45 GLT28D1

Immunogen:

Recombinant Human Putative bifunctional UDP-N-acetylglucosamine transferase and deubiquitinase ALG13 protein (112-285AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A