

# DIAPH3 Antibody, FITC conjugated



PACO49328

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds to GTP-bound form of Rho and to profilin. Acts in a Rho-dependent manner to recruit profilin to the membrane, where it promotes actin polymerization. It is required for cytokinesis, stress fiber formation, and transcriptional activation of the serum response factor. DFR proteins couple Rho and Src tyrosine kinase during signaling and the regulation of actin dynamics.

**Gene ID:**

DIAPH3

**Uniprot**

Q9NSV4

**Synonyms:**

Protein diaphanous homolog 3 (Diaphanous-related formin-3) (DRF3) (MDia2), DIAPH3, DIAP3

**Immunogen:**

Recombinant Human Protein diaphanous homolog 3 protein (642-849AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A