

# USE1 Antibody, HRP conjugated



PACO49343

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.

**Gene ID:**

USE1

**Uniprot**

Q9NZ43

**Synonyms:**

Vesicle transport protein USE1 (Putative MAPK-activating protein PM26) (USE1-like protein) (p31), USE1, USE1L

**Immunogen:**

Recombinant Human Vesicle transport protein USE1 protein (1-146AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A