

PCMP-H63 Antibody



PACO49350

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PCMP-H63

Uniprot

Q9SUU7

Synonyms:

Pentatricopeptide repeat-containing protein At4g32450, mitochondrial, PCMP-H63

Immunogen:

Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At4g32450, mitochondrial protein (111-537AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A