

AMPD2 Antibody, HRP conjugated



PACO49395

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

AMP deaminase plays a critical role in energy metabolism. Catalyzes the deamination of AMP to IMP and plays an important role in the purine nucleotide cycle.

Gene ID:

AMPD2

Uniprot

Q01433

Synonyms:

AMP deaminase 2 (EC 3.5.4.6) (AMP deaminase isoform L), AMPD2

Immunogen:

Recombinant Human AMP deaminase 2 protein (1-129AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A