

# AMPD2 Antibody, FITC conjugated



PACO49396

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

AMP deaminase plays a critical role in energy metabolism. Catalyzes the deamination of AMP to IMP and plays an important role in the purine nucleotide cycle.

**Gene ID:**

AMPD2

**Uniprot**

Q01433

**Synonyms:**

AMP deaminase 2 (EC 3.5.4.6) (AMP deaminase isoform L), AMPD2

**Immunogen:**

Recombinant Human AMP deaminase 2 protein (1-129AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A