

CALCA Antibody, FITC conjugated



PACO49444

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Causes a rapid but short-lived drop in the level of calcium and phosphate in blood by promoting the incorporation of those ions in the bones.

Gene ID:

CALCA

Uniprot

P41547

Synonyms:

Calcitonin, CALCA, CALC CCALCI

Immunogen:

Recombinant Dog Calcitonin protein (85-116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A