

# EN1 Antibody



PACO49562

---

## Product Information

**Size:**

50ug

**Reactivity:**

Chicken

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EN1

**Uniprot**

Q05916

**Synonyms:**

Homeobox protein engrailed-1 (Gg-En-1) (Homeobox protein en-1), EN1, EN-1

**Immunogen:**

Recombinant Chicken Homeobox protein engrailed-1 protein (1-333AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A