

EN2 Antibody, Biotin conjugated



PACO49569

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleus, multicellular organism development.

Gene ID:

EN2

Uniprot

P19622

Synonyms:

Homeobox protein engrailed-2 (Homeobox protein en-2) (Hu-En-2), EN2

Immunogen:

Recombinant Human Homeobox protein engrailed-2 protein (1-251AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A