

# GPR85 Antibody, HRP conjugated



PACO49651

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Orphan receptor.

**Gene ID:**

GPR85

**Uniprot**

P60893

**Synonyms:**

Probable G-protein coupled receptor 85 (Super conserved receptor expressed in brain 2), GPR85, SREB2

**Immunogen:**

Recombinant Human Probable G-protein coupled receptor 85 protein (229-267AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A