

HMBS Antibody, HRP conjugated



PACO49655

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Tetrapolymerization of the monopyrrole PBG into the hydroxymethylbilane pre-uroporphyrinogen in several discrete steps.

Gene ID:

HMBS

Uniprot

P08397

Synonyms:

Porphobilinogen deaminase (PBG-D) (EC 2.5.1.61) (Hydroxymethylbilane synthase) (HMBS) (Pre-uroporphyrinogen synthase), HMBS, PBGD UPS

Immunogen:

Recombinant Human Porphobilinogen deaminase protein (283-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A