

HOXD12 Antibody



PACO49662

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXD12

Uniprot

P35452

Synonyms:

Homeobox protein Hox-D12 (Homeobox protein Hox-4H), HOXD12, HOX4H


Immunogen:

Recombinant Human Homeobox protein Hox-D12 protein (30-169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO49662 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).