

HOXD12 Antibody, HRP conjugated



PACO49663

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXD12

Uniprot

P35452

Synonyms:

Homeobox protein Hox-D12 (Homeobox protein Hox-4H), HOXD12, HOX4H

Immunogen:

Recombinant Human Homeobox protein Hox-D12 protein (30-169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A