

# HOXD12 Antibody, Biotin conjugated



PACO49665

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

**Gene ID:**

HOXD12

**Uniprot**

P35452

**Synonyms:**

Homeobox protein Hox-D12 (Homeobox protein Hox-4H), HOXD12, HOX4H

**Immunogen:**

Recombinant Human Homeobox protein Hox-D12 protein (30-169AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A