

PACO49702

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to one of the cis-acting domain of the insulin gene enhancer. May be involved in subtype specialization of primary motoneurons.

Gene ID:

isl1

Uniprot

P53405

Synonyms:

Insulin gene enhancer protein isl-1 (Islet-1), isl1, isl-1

Immunogen:

Recombinant Zebrafish Insulin gene enhancer protein isl-1 protein (1-349AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A