

MLPH Antibody, FITC conjugated



PACO49776

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Rab effector protein involved in melanosome transport. Serves as link between melanosome-bound RAB27A and the motor protein MYO5A.

Gene ID:

MLPH

Uniprot

Q9BV36

Synonyms:

Melanophilin (Exophilin-3) (Slp homolog lacking C2 domains a) (Slac2-a) (Synaptotagmin-like protein 2a), MLPH, SLAC2A

Immunogen:

Recombinant Human Melanophilin protein (310-539AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A